

# Anaplasma phagocytophilum

## PERFORMANCE CHARACTERISTICS

### Sensitivity

The indirect immunofluorescence antibody assay (IFA) for *Anaplasma phagocytophilum* has been used for over 25 years. In the case of naturally-infected animals and humans, the appearance of diagnostic antibody often lags behind the appearance of clinical signs, with the majority of (acute) sera negative at presentation. Within 1-2 weeks most subsequently seroconvert to high titers. The first Fuller Laboratories IFA product was introduced in 1995 and utilized the MRK equine strain grown in KG-1 cells. Comparison between this former substrate and the current NCH-1 isolate grown in HL60 cells (25 positive and 25 negative sera) demonstrates 100% concordance. Correlations of IFA protocols with Western Immunoblot (WB) techniques demonstrate IFA sensitivity between 80-100% <sup>3-</sup>. An inhouse series of 48 dog sera from the state of New York showed complete (100%) concordance between IFA and WB on 20 positive sera and 28 negative.

### Specificity

With the incorporation of *Ehrlichia equi*, *Ehrlichia phagocytophilum* and the HE Agent into the species *Anaplasma phagocytophilum* comb. nov., there are few close relatives to this combined species. By IFA there have been reports of human serum crossreactivity with *Ehrlichia chaffeensis*. There are no sources of crossreactivity outside the tribe Ehrlichiae. Canine sera from a non-endemic region, metropolitan Southern California, were tested as a source of negative sera. Of 58 sera tested, all (58/58) were negative (100% specificity).

## REFERENCES

1. Greig, B., Asanovich, K.M., Armstrong, P.J. and Dumler, J.S. 1996. Geographical, Clinical, Serologic, and Molecular Evidence of Granulocytic Ehrlichiosis, a Likely Zoonotic Disease, in Minnesota and Wisconsin Dogs. *J. Clin. Microbiol.* 34:44- 48.
2. Massung, R.F., K. Slater, J.H. Owens, W.L. Nicholson, T.M. Mather, V.B. Solderg, and J.G. Olson. 1998. Nested PCR Assay for Detection of granulocytic Ehrlichiae. *J. Clin. Microbiol.* 36:10901095.
3. Zhu, N., N. Ohashi, Y. Rikihisa, H.W. Horowitz, G.P. Wormser, and K. Hechemy. 1998. Cloning and Expression of the 44- Kilodalton Major Outer Membrane Protein Gene of the Human Granulocytic Ehrlichiosis Agent and Application of the recombinant Protein to Serodiagnosis. *J. Clin. Microbiol.* 36:1666-1673.
4. Zhu, N., Y. Rikihisa, H.Y. Kim, G.P. Wormser, and H.W. Horowitz. 1997. Comparison of Major Antigenic Proteins of Six Strains of the Human Granulocytic Ehrlichiosis Agent by Western Immunoblot Analysis. *J. Clin. Microbiol.* 35: 2606-2611.

5. Breitschwerdt, E.B., Hegarty, B.C. and Hancock, S.I. 1998. Sequential Evaluation of Dogs Naturally Infected with Ehrlichia canis, Ehrlichia chaffeensis, Ehrlichia equi, Ehrlichia ewingii, or Bartonella vinsonii. J. Clin. Microbiol. 36:265-2651.

GENTAUR